

Baxter **HD**

Now from Baxter...  
the EXELTRA family of  
high-flux, single-use dialyzers

**extraordinary clearance**  
**ultra compatibility**



**EXELTRA**  
SINGLE-USE DIALYZERS

**Baxter**

# challenge

other dialyzers with these performance comparisons



## EXELTRA Dialyzer Clearance Data

	Urea				Creatinine				Phosphate				Vitamin B <sub>12</sub>			
	Qb (mL/min)	200	300	400	500	200	300	400	500	200	300	400	500	200	300	400
Qd (mL/min)	500				500				500				500			
<b>EXELTRA 150</b>	193	262	305	332	186	242	274	297	179	227	255	274	132	152	163	170
<b>EXELTRA 170</b>	196	268	310	341	190	252	286	307	179	232	261	280	138	160	172	180
<b>EXELTRA 190</b>	197	273	323	354	190	251	289	313	186	242	276	296	143	168	183	193
<b>EXELTRA Plus 210</b>	199	287	350	384	198	277	328	363	191	252	292	318	164	202	222	232

## EXELTRA Dialyzer Technical Information

<b>Membrane Polymer</b>	Cellulose triacetate
<b>Sterilization</b>	Gamma irradiation
<b>Membrane Inner Diameter</b>	200 µm
<b>Membrane Thickness</b>	15 µm
<b>Maximum TMP</b>	500 mm Hg (66 kPa)

## EXELTRA Dialyzer Specifications

	EXELTRA 150	EXELTRA 170	EXELTRA 190	EXELTRA Plus 210
<b>Code Number</b>	5M2119	5M2120	5M2121	5M2132
<b>Surface Area (m<sup>2</sup>)</b>	1.5	1.7	1.9	2.1
<b>Priming Volume (mL, aqueous)</b>	95	105	115	125
<b>Ultrafiltration Rate (mL/hr/100 mm Hg)</b>	3150	3380	3642	4736
<b>KoA (Qb 300, Qd 500)</b>	993	1103	1214	1714

To find out how the Baxter EXELTRA family of dialyzers can benefit your patients, contact your Baxter Account Executive at 1-888-736-2543, option 4.

**EXELTRA**  
 SINGLE-USE DIALYZERS  
*Extraordinary Clearance.  
 Ultra Compatibility.*

**Baxter**

Baxter Healthcare Corporation  
 Renal Division  
 1620 Waukegan Road  
 McGaw Park, IL 60085  
 1-888-736-2543

Baxter and EXELTRA are trademarks  
 of Baxter International Inc.

©Copyright 2003, Baxter Healthcare Corporation  
 5L0467 10/03  
 AL03374