



Hemodialysis Concentrates

Baxter and Rockwell Medical

Hemodialysis Concentrates

Baxter Expands Portfolio: Baxter and Rockwell Medical

Baxter is now the exclusive distributor of Rockwell's hemodialysis concentrate products in the US. Rockwell is a leading manufacturer of products targeting chronic kidney disease with innovative products and services. This collaboration expands Baxter's product portfolio with the addition of Rockwell's hemodialysis concentrate products.



Citric acid concentrate powders and liquids

CitraPure® acid concentrate is formulated to be used in conjunction with RenalPure® Powder Bicarbonate (CB-Series) or SteriLyte® Liquid Bicarbonate (CL-Series) only in a 45X three-stream hemodialysis machine. The listed electrolyte values represent the final dialysate electrolyte yield.

Formula	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Citric Acid (mEq/L)	Dextrose (mg%)
45X CitraPure Liquid Acid Concentrate – Drum (55-gallon)						
L4-108-55	137	2.00	2.00	1.00	2.40	100
L4-116-55	137	2.00	2.50	1.00	2.40	100
L4-117-55	137	3.00	2.50	1.00	2.40	100
45X CitraPure Liquid Acid Concentrate – Case (four one-gallon containers/case)						
L4-108-4	137	2.00	2.00	1.00	2.40	100
L4-116-4	137	2.00	2.50	1.00	2.40	100
L4-117-4	137	3.00	2.50	1.00	2.40	100
45X CitraPure Dry Acid Concentrate – Case (each case produces 25 gallons of CitraPure Liquid)						
D4-108	137	2.00	2.00	1.00	2.40	100
D4-116	137	2.00	2.50	1.00	2.40	100
D4-117	137	3.00	2.50	1.00	2.40	100

CitraPure®

Acetic acid concentrate powders and liquids

RenalPure® and Dri-Sate® acid concentrates are formulated to be used in conjunction with RenalPure Powder Bicarbonate (CB-Series) or SteriLyte® Liquid Bicarbonate (CL-Series) only in a 45X three-stream hemodialysis machine. The listed electrolyte values represent the final dialysate electrolyte yield.

Formula	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Cl- (mEq/L)	Acetic Acid (mEq/L)	Dextrose (mg%)
45X RenalPure Liquid Acid Concentrate – Drum (55-gallon)							
C-108-55	137	2.00	2.00	0.75	105	4.00	100
C-116-55	137	2.00	2.50	0.75	105	4.00	100
C-117-55	137	3.00	2.50	0.75	106	4.00	100
45X RenalPure Liquid Acid Concentrate – Case (four one-gallon containers/case)							
C-108-4	137	2.00	2.00	0.75	105	4.00	100
C-116-4	137	2.00	2.50	0.75	105	4.00	100
C-117-4	137	3.00	2.50	0.75	106	4.00	100
45X Dri-Sate Dry Acid Concentrate – Case (each case produces 25 gallons of liquid acid)							
DC-108	137	2.00	2.00	0.75	105	4.00	100
DC-116	137	2.00	2.50	0.75	105	4.00	100
DC-117	137	3.00	2.50	0.75	106	4.00	100

Bicarbonate concentrate powders and liquids

RenalPure Powder Bicarbonate and SteriLyte Liquid Bicarbonate are for use with CitraPure® Dry or Liquid Acid Concentrate, Dri-Sate Dry Acid or RenalPure Liquid Acid Concentrate.

Formula	Na+ (mEq/L)	Bicarbonate (mEq/L)	Qty/case	Volume produced per bag (gallons)
RenalPure Powder Bicarbonate Bags				
CB-02-20C	37.0	37.0	20	2.1
CB-15-4C	37.0	37.0	4	15.0
CB-25-2C	37.0	37.0	2	25.0
SteriLyte Liquid Bicarbonate – Case (four one-gallon containers/case)				
CL-02-4	37.0	37.0	4	NA

Custom formulations for 45X are available upon request. Other formulations for 35X and 36.83X are also available upon request. Contact your sales representative for more information.

CitraPure acid concentrate, Dri-Sate acid concentrate, RenalPure acid concentrate, RenalPure Powder Bicarbonate and SteriLyte Liquid Bicarbonate are indicated for use in acute

and chronic hemodialysis and to be used with the appropriate three stream hemodialysis machine/monitor and bicarbonate concentrate.

Rx Only. For the safe and proper use of the concentrates mentioned herein, refer to the *Instructions for Use*.

Contact your local Baxter sales representative to learn more about the product.
For Customer Support call 800-525-2623.

Baxter is a trademark of Baxter International Inc.
Citrapure, Dri-Sate, RenalPure, and SteriLyte are trademarks of Rockwell Medical.
www.baxter.com